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
PELVIPOWER™

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MAGNETFELD-TRAINING

PELVIPOWER ACADEMY

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
Voiding Dysfunction

Pulsed Magnetic Stimulation for Stress Urinary Incontinence: 1-Year Followup Results

Renly Lim,^{*,†} Men Long Liang, Wing Seng Leong, Nurzalina Abdul Karim Khan and Kah Hay Yuen

From the School of Pharmaceutical Sciences, Universiti Sains Malaysia (USM, NAKK, KHI) and Departments of Urology, Island Hospital (MLI) and Lam Wah Ee Hospital (WSE), Penang, Malaysia, and School of Pharmacy and Medical Sciences, University of South Australia (USA), Adelaide, South Australia, Australia

Abbreviations and Acronyms
ICI = International Consultation on Incontinence
ICIQ-LUTS(q) = ICI Questionnaire-Lower Urinary Tract Symptoms Quality of Life
ICIQ-UI SF = ICI Questionnaire for Urinary Incontinence-Short Form
PFM = pelvic floor muscle
PFMT = PFM training
PGI-I = Patient Global Impression of Improvement
PMS = pulsed magnetic stimulation
SUI = stress urinary incontinence

Purpose: Despite significant differences in success rates between surgical and nonsurgical treatments for female stress urinary incontinence, a few cross-sectional surveys showed that most patients still prefer the latter. We evaluated the efficacy of the under studied nonsurgical treatment using pulsed magnetic stimulation for female stress urinary incontinence.

Materials and Methods: This randomized, double-blind, sham controlled study was performed in 120 female subjects at least 21 years old with stress urinary incontinence. Treatment involved pulsed magnetic stimulation for 2 sessions per week for 2 months (16 sessions). After 2 months, subjects could opt for 16 additional sessions regardless of initial randomization. The primary response criterion was a 5-point reduction in the ICIq-UI SF (International Consultation on Incontinence Questionnaire for Urinary Incontinence-Short Form) score. Key secondary response criteria included objective and subjective cure, supplemented by other secondary criteria. Followups were performed at months 1, 2, 5, 8 and 14.

Results: At 2 months 45 of 60 subjects (75%) in the active arm vs 13 of 60 (21.7%) in the sham arm were treatment responders ($p < 0.001$). After 2 months 24 subjects (40%) in the active arm and 41 (68%) in the sham arm elected additional active pulsed magnetic stimulation. At 14 months, subjects who received 32 sessions of active pulsed magnetic stimulation had the highest percentage of treatment responders (18 of 24 or 75.0%), followed by those who received 16 sessions (26 of 36 or 72.2% and 28 of 41 or 68.3%) and those who did not receive any active pulsed magnetic stimulation (4 of 19 or 21.1%) ($p < 0.001$).

Conclusions: The encouraging long-term response rates show that pulsed magnetic stimulation is an attractive nonsurgical alternative for patients who do not want to undergo surgery.

Key Words: urinary bladder; urinary incontinence, stress; pelvic floor; magnetic field therapy; risk

STRESS urinary incontinence is a common and distressing condition.^{1,2} The 5th ICI advocated PFMT with success rates of 15% to 56% as the gold standard nonsurgical treatment.^{3,4} However, there is no

standardized PFMT regimen⁵ and the success of PFMT is often limited by poor compliance.^{6,7} Other nonsurgical options (eg biofeedback, vaginal cones and electrical stimulation) are limited by low success rates of 9% to 65%,

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